

Fully automated microscope incubation solutions for high content live cell imaging applications



Complete control
never stops, even in motion.


BOLD **LINE 3**





Motorized door/lid

Externally triggered doors allow access to robotic plate handler for stage top or cage configurations.



Precision repositioning

Motorized plunger ensures precise plate repositioning every time, and loading sensor provides status for robotic coordination.




Plate hotel fully incubated

Compatible with robotic loader, this enclosure provides heat, CO₂ and Humidity to store multiple plates during automated microscopy.



Auto refill for humidifier

Never worry about running out of water; refill maintains digital humidity control for weeks without user intervention.



CUSTOM SOLUTIONS AVAILABLE

At Okolab, we take pride in our ability to create custom solutions uniquely suited for our biotech customers' most demanding applications. Any of the above solutions can be custom-made or modified to best fit your automated workflow.

Have an incubation need for your automated system that isn't listed here?

Whether it is for custom-designed, heated enclosures for non-standard setups or for specialized gas control requirements, **let us design and build it for you.**

oko-lab.com

