

The **ONLY** incubator with total sample feedback





-  CO₂
-  Humidity
-  Temperature

Complete control
where it matters most.

BOLD **LINE 3**

In-chamber feedback



Total Sample Control

In-Chamber and closed-loop CO₂, Heating, and Humidity feedback, ensuring improved accuracy for precise and reliable experiments.



Ease of Use

Improved user experience with an optimized interface and an efficient ductless design that minimizes footprint in your laboratory space.




Continued Quality

Based on 20 years of incubation experience and backed by an exclusive 3-year warranty on core components.



Fast and Stable

Highly responsive, double heaters provide unprecedented temperature uniformity with a 50% faster warmup time.



NEW UNIVERSAL INSERTS

Convenient design allows our sample inserts to be used outside of our chambers as standard, universal sample holders for fixed sample imaging on virtually any microscope stage.

oko-lab.com

